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NANOREPRO SARS-Cov-2 Antigen Rapid Test

INTENDED USE

This rapid test kit is intended for the qualitative detection of SARS-Cov-2 infection from patients. It is for professional use only. It is an aid in the diagnosis of the patients with suspected SARS-Cov-2 infection in conjunction with clinical presentation and results of other laboratory tests. Results from this test kit should not be used as the sole basis for diagnosis.

The test provides preliminary test results. Negative results do not preclude SARS-Cov-2 infection and should not be used as the sole basis for treatment or other management decision.

PRINCIPLE

This kit is an immunochromatography assay. According to the gold immunochromatographic test principle, double antibody sandwich method was used to detect SARS-Cov-2 antigen in the samples. When there is virus antigen present in the sample, the antigen binds with the corresponding colloidal gold monoclonal antibody and the coated monoclonal antibody at the detection line to form a compound and then condenses into a red band, indicating a positive result. If there is no antigen in the sample, complex cannot be formed at the detection line, and no red band is shown, indicating negative result.

Whether the sample contains antigen or not, the gold monoclonal antibody will bind to the enveloped antibody at the quality control line, form a compound and condense into a red band.

KIT COMPONENTS

- 20 Test cassettes
- 2 Sample extraction
- 20 Sample tubes
- 20 single packaged Swabs
- 1 Tube stand
- 1 Instruction for use

STORAGE AND EXPIRY

- Store as packaged in the sealed pouch at 2-30°C, avoid hot and sunshine, dry place, valid for 12 months.
- Do not freeze. Some protective measures should be taken in hot summer and cold winter to avoid high temperature or freeze-thaw.
- Do not open the inner packaging until ready, it must be used in one hour if opened (Humidity≤60%, Temp: 20°C-30°C).
- Please use immediately when the humidity>60%.

SAMPLE PRESERVATION

Sample of human nasopharyngeal swabs and oropharyngeal swabs should be processed as soon as possible after sample collection. If the test cannot be performed immediately, the sample should be stored in a sealed state, stored at 2-8°C for 8 hours, and stored below -20°C for 1 month. Long-term storage is not recommended.

SAMPLE COLLECTION

By nasopharyngeal swab: The sampler gently holds the head of the person to be collected with one hand; holds the swab with the other hand, sticks the swab to the nostril to enter, and slowly penetrates backwards along the bottom of the lower nasal passage, so as not to exert too much force to avoid traumatic hemorrhage. When the tip of the swab reaches the posterior wall of the nasopharynx cavity, gently rotate it (make sure the swab contains both cells and slime), and then slowly take out the swab.


By oropharyngeal swab: The head of the person to be collected is slightly tilted and his mouth is wide open, exposing the pharyngeal tonsils on both sides. Wipe the swab across the root of the tongue.

Wipe the pharyngeal tonsils on both sides of the person to be collected back and forth with a little force for at least 3 times, and then wipe up and down the posterior pharyngeal wall for at least 3 times.


TEST STEPS

Instructions must be read entirely before taking the test. Bring test cassettes and test components at room temperature before performing the test.

1. Add 23 drops sample extraction (550 µl) into the sample tube.

2. Dip the swab after collecting the sample into the sample extraction liquid, fully immerse the tip of the swab, rotate and squeeze the swab 10 times. Roll the swab head against the inner wall of the tube as you remove it. Try to release as much liquid as possible. Dispose the used swab in accordance with biohazard waste disposal protocol.

3. Take out the test cassette from the sealed pouch, place it on a clean and level surface with the sample port well up. Name the test cassette.

5. Observe the test results immediately within 15~20 minutes, the result is invalid over 20 minutes.

Note: When using viral transport medium (VTM), it is important to ensure that the VTM containing the sample is warmed to room temperature. Cold samples will not flow correctly and can lead to erroneous or invalid results. Several minutes will be required to bring a cold sample to room temperature.

INTERPRETATION OF RESULTS

POSITIVE:

Two (2) distinct colored lines appear. One line should be in the control region (C) and the other line should be in the test region (T).


NEGATIVE:

One (1) colored line appears in the control region (C). No apparent colored line appears in the test region (T). The negative result does not indicate the absence of analytes in the sample.

INVALID:

No colored lines appear, or control line fails to appear, indicating that the operator error or reagent failure. Verify the test procedure and repeat the test.

2. Limit of detection:

LOD concentration	30 TCID ₅₀ /mL
-------------------	---------------------------

3. Cross-reactivity

With human coronavirus 229E, human coronavirus OC43, human coronavirus NL63, adenovirus, human metapneumovirus, parainfluenza virus 1, parainfluenza virus 2, parainfluenza virus 3, parainfluenza virus 4, influenza A, type B influenza, enterovirus, respiratory syncytial virus, rhinovirus, SARS coronavirus, MERS coronavirus, Haemophilus influenzae, Streptococcus pneumoniae, Streptococcus pyogenes, Candida albicans, Bordetella pertussis, Mycoplasma pneumoniae, Legionella pneumophila, etc. have no cross reaction (see table in the German IFU).

4. Interfering

Common interfering substances in the sample, such as blood, mucin, and pus, have no effect on the test results.

PRECAUTIONS

1. For IN VITRO diagnostic use only.
2. Reagents should be used as soon as possible after opened. This reagent cannot be reused for disposable.
3. The test device should remain in the sealed pouches until use.
4. All specimens and reagents should be considered potentially hazardous and handled in the same manner as an infectious agent after use.
5. Don't use after the expiration date.

LITERATURHINWEISE | REFERENCES

1. Formi, D., Cagliani, R., Clerici, M. & Sironi, M. Molecular evolution of human coronavirus genomes. Trends Microbiol. 25, 35-48 (2017).
2. Ithete, N. L. et al. Close relative of human Middle East respiratory syndrome coronavirus in bat, South Africa. Emerg. Infect. Dis. 19, 1697-1699 (2013).

PERFORMANCE CHARACTERISTICS

1. **Sensitivity and Specificity**

The verification of naso-oropharyngeal samples is carried out with enterprise reference and the results meet the following requirements:

The National Reference Panel for SARS-Cov-2 Antigen Detection Kit of the National Institutes for Food and Drug Control was used for verification, and the results met the following requirements:

Coincidence rate of negative reference products: National negative reference materials (N1-N20) were tested, and the results were all negative, with a coincidence rate of 100%. Coincidence rate of positive reference products: tested with national positive reference products (P1-P8), the results are all positive, and the coincidence rate is 100%.

Repeatability: The 10 test results of R1 and R2 should be positive, and the color rendering is uniform.

Clinical Sensitivity and Specificity:

SARS-Cov-2 Antigen Rapid Test	RT-PCR		Total
	Positive	Negative	
Positive	243	9	252
Negative	7	541	548
Total	250	550	800

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NanoRepro AG
Nanorepro AG