



Trusted by families for over 40 years

# User Manual

## NON-MEDICAL USE ISOLATION GOWN (LVL-3)

One Size Fits All (Large / X Large)

### Class 1 Device

#### Product Specification & Size:

Unit: Inches

Size	Length(From the shoulder to bottom)	Width(Single Layer)	Belt	Sleeve Length	Tolerance
(S) 42.5" x 54.3"	43.3"	25.6"	24"	24"	±0.8"
(M) 45.3" x 55.1"	47.2"	27.6"	24"	24.4"	±0.8"
(L) 47.2" x 57.1"	51.2"	27.6"	24"	24.4"	±0.8"
(XL) 49.2" x 59.1"	55.1"	29.5"	24"	24.8"	±0.8"
(XXL) 54.3" x 63"	59.1"	31.5"	24"	24.8"	±0.8"


**Product Structure:** The product is made up of non-woven fabric material with certain protective properties and is composed of front body, back body, sleeves and ties.

**Fabric Weight:** 55 GSM

#### Material Content:

PP + PE (55% polypropylene + 45% polyethylene)

**Fabric Composition:** 30GSM SBPP non woven laminated with 25GSM PE film. Extra strength double layer protection.



Front

Back

#### Product Intended Use

- This disposable isolation gown is intended to protect the wearer from the transfer of microorganisms and bodily fluids in low-or-minimal risk patient isolation situations.

#### Features:

- Durable, breathable, and odorless
- Comfortable fit
- Highly resistant and protective against water and fluids
- Latex free
- All seams have inner binding
- Ties around neck and waist for extra protection

**FOR REORDERS CONTACT:** [sales@homewearlinens.com](mailto:sales@homewearlinens.com)

Distributed by Sam Hedaya Corporation  
295 Fifth Ave, New York, NY 10016

### Storage and Disposal

- Store in dry, clean conditions in original packaging, away from direct sunlight sources of high temperature and solvent vapors.
- Store within the temperature range of -30°C to +55°C (-22°F to +131°F) and with relative humidity below 80%.
- Expected shelf life is three years from date of manufacture when stored as directed.
- Replace garments if damaged, heavily contaminated or in accordance with applicable regulations.

### Limited Use



Do not wash



Do not dry clean



Do not bleach



Flammable - keep away from sparks or flames



Do not iron



Single Use - do not re-use



Do not tumble dry

### Warnings and Limitations

Before use, read and understand all user instructions and be sure that the product is suitable for the application and fits correctly. Product must never be altered or modified.

#### Do not use for:

- Contact with heavy oils, sparks, or flame or combustible liquids
- Environments with high mechanical risks (abrasive surfaces, sharp edges, etc.)
- Contact with hazardous substances beyond type 4/5/6 certification
- Environments with conditions of excessive heat
- Environments with a presence of flammable gas
- Any surgical setting or where significant exposure to liquid bodily or other hazardous fluids may be expected

Manufacturer: Zhejiang Huafu Medical Equipment Co., Ltd  
ADDRESS: 688 Xinxing 1st RD., Pinghu, Jiaxing  
Zhejiang Province 314200 China

40 mm



Trusted by families  
for over 40 years

**PP/PE ISOLATION GOWN**  
**NON-MEDICAL USE**  
**ONE SIZE FITS ALL/X-LARGE**  
**10pcs/bag**

400 mm

320 mm

**WARNING:**  
**TO AVOID DANGER OF SUFFOCATION, KEEP THIS**  
**BAG AWAY FROM BABIES AND CHILDREN. DO NOT USE IN CRIBS, BEDS,**  
**CARRIAGES OR PLAYPENS. THIS BAG IS NOT A TOY.**

**AVERTISSEMENT:**  
**POUR EVITER TOUT DANGER DE**  
**SUFFOCATION, GARDEZ CE SAC HORS DE LA PORTÉE DES BÉBÉS ET DES ENFANTS.**  
**NE PAS UTILISER DANS UN LIT D'ENFANT, LITS LANDAUS**  
**OU PARCS POUR BÉBÉS, CE SAC N'EST PAS UN JOUET**

**PRECAUCION:**  
**PARA EVITAR EL PELIGRO DE ASFIXIA,**  
**MANTENGA ESTA BOLSA LEJOS DE BEBÉS Y NIÑOS. NO LA USE EN**  
**CUNAS, CAMAS, CARREOLAS O CORRALES DE JUEGO, ESTA BOLSA**  
**NO ES UN JUGUETE.**

"GOWN" - PP/PE  
Medical Isolation Gown/Disposable



[sales@homewearlinens.com](mailto:sales@homewearlinens.com)

400 mm

320 mm